CareMat

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Designed in Switzerland / Made in Switzerland

CareMat

with 434 MHz TeleAlarm transmitter

The pressure-sensitive bed mat
for dementia patients and people at risk of falling

Translation of the original instructions

1 Safety instructions

- Read these operating instructions thoroughly before putting the device into operation and keep them for future reference.
- Not a safety component in accordance with the EU Machinery Directive must not be used for personal protection or EMERGENCY STOP function.
- Do not use this product other than for its specified application.
- Failure to follow these safety precautions may cause damage to device or objects, serious personal injury, or death.
- Observe all applicable local, national, and international safety standards, codes, and laws.
- After accessing the inside of the device, ensure the cover/protection seal is closed tightly to achieve designated protection rating.
- Warranty claims will be invalid in the event of any damage to the product, to goods or persons that is due to failure to observe these instructions.

Instructions for nursing staff

- The CareMat operates purely as a support for the nursing staff, but is not a substitute for them.
- The functionality of the CareMat is to be checked each time that it is used.

Notes for the installer

- The electrical installation work and commissioning of the product may only be carried out by a qualified technician.
- Avoid touching any electronic components.
- ESD safety rules should be observed.

2 Product description


2.1 Described radio system

The CareMat with 434 MHz Radio Transmitter is compatible with the LE10 Radio Receiver from TeleAlarm. The LE10 can receive radio signals from different TeleAlarm radio transmitters and relay them to a cable connected call system. Detailed information on this and further compatible 434 MHz appliances is available on the following website:
<http://www.telealarm.com/en/central-components>


2.2 Areas of application and function

The CareMat was primarily designed for use in elderly people's and care homes, hospitals and psychiatric wards as well as for sheltered residential accommodation. The necessary care and supervision of patients with dementia and of those who are at risk of falling represents a significant increase in care expenditure for care practitioners. The CareMat is the ideal support: Risks of falling and the dangers which arise from wandering or loss of orientation can be significantly reduced.


The CareMat, as a pressure-sensitive sensing device, is laid directly in front of the care bed. One LE10 transmitter per room is to be installed and connected to the internal call system. As soon as the person leaves the bed and steps onto the CareMat, it switches on and sends out an alarm signal to the call system via the LE10. The wireless signalling system reduces the hazardous tangle of wiring around the care bed.

CareMat A01T B434
CareMat A11T B434
1100 x 700 mm

CareMat B01T B434
CareMat B11T B434
700 x 400 mm


CareMat C01T B434
CareMat C11T B434
1100 x 700 mm,
semi-circular

- ① Housing with 434 MHz transmitter
- ② Chamfered, wheelchair-accessible edges
- ③ Non-slip surface with stud texture design


- ④ Carrying handle
- ⑤ Active switching zone
- ⑥ 434 MHz LE10 signal receiver (optional component)

3 Day to day operation

The CareMat is centred and placed in front of the bed with the handles underneath the edge of the bed to ensure that missing the mat is impossible.

Tune-in the CareMat to the LE10 signal receiver when using for the first time.

Should the CareMat not be needed, it can be pushed underneath the bed.

4 Cleaning and disinfecting

- The CareMat can be cleaned using commercially available disinfecting and cleaning agents with max. 70 vol% alcohol. No corrosive agents should be used.
- The CareMat is neither autoclavable nor can it be disinfected using steam.

5 Transport and storage

- During transport and storage, any significant bending of the CareMat should be avoided.
- Store the CareMat in a dry, clean and flat location.
- Do not place any heavy loads on the CareMat and do not stack more than 10 packed mats on top of each other.


6 Signal system combined with light call system

6.1 Configuration LE10-radio receiver


One LE10 transmitter, which is connected to the call system via the room terminal, is to be installed per room. The factory settings are to be retained: all S3 micro switches are to be OFF. This corresponds to the Normal operating mode. The LED 1 (green) is constantly lit in the Standby mode. If any-one steps on an assigned CareMat, an alarm signal is transmitted to the LE10. The LED 2 (red) illuminates and the Relay 1 output switches on for 10 seconds.

A run-down transmitter battery (charge condition < 20%) is displayed in the following manner: LED 2 (red) flashes every 10 seconds and Relay 2 is switched on. When the battery has been changed, LED 2 and Relay 2 are switched-off as soon as a daily data transmission or an activation of the CareMat has been recognised.

LE10-radio receiver

Connector compartment

Inputs

The LE10 can be supplied in two ways:

A) via the connections 1 and 2 of the P130 cable connector or

B) via the RJ12-jack for plug-in power supply use < 25 mA

Output Relays

max. 28 V AC/DC

max. 250 mA

6.2 Programming the 434 MHz radio transmitter

Tuning the CareMat to the LE10: Press "T1" key for 2 seconds. LED 2 (red) flashes slowly. Tune-in CareMat within 30 seconds by activating it. When the programming has been successfully completed, LED 1 and LED 2 flash alternately. If the CareMat has already been programmed, LED 2 flashes rapidly for six seconds. When not being used or when being used in another room, the CareMat should be switched off on the LE10. This way it can be ensured that the call is only triggered in the specific room where the CareMat is actually located.


To switch off CareMat on LE10: Press the "T1" key for 5 seconds until the LED 1 (green) starts to flash. Then activate CareMat and wait for 5 seconds. The CareMat is now switched-off.

Note: If no CareMat, or no transmitter is tuned-in on the LE 10, the LED 2 (red) flashes three times every 10 seconds. In order to avoid this, a CareMat should always be tuned-in on the LE10.

6.3 Battery changing

Remove housing cover on rear
(6 countersunk screws M3x5, cross-head size PH1)

- ① With a thin, pointed object (preferably plastic), push the battery out of the housing compartment.
- ② Insert the new battery, ensuring the positive (+) side is facing you (as shown in the drawing above).

Reinstall housing cover, carefully tighten the 6 screws by hand, do not overtighten.

7 Maintenance and fault finding

The CareMat is maintenance-free up until the time that the batteries need changing. Only tighten the screws lightly!

If despite this the function is not being provided, the following checks should be made:

- Is the call system operating and the LE10 connected correctly?
- Is the CareMat correctly tuned-in to the LE10?
- Is the signal transmitter battery OK?

If the points mentioned above have been checked out, but the functions of the CareMat are, despite this, still not being provided, you should consult a qualified expert.

8 Technical specifications

Article	CareMat Ax1T B434 / Bx1T B434 / Cx1T B434	Signaling protocol	TeleAlarm
Dimensions (LxWxH)	1100x700x9 mm / 700x400x9 mm / 1100x700x9 mm, height with housing 15 mm	Radio frequency	434.01 MHz
Weight [kg]	8.2 kg / 3.2 kg / 6.5 kg	Working life	3 years (with 3 alarms and 1 test-status per day)
Material	Polyurethane (PU)	Radio range	up to 30 m in buildings, dependent on structural building conditions (up to 250 m in free field)
Colour	Careyellow, Caregrey	Type of battery	1x CR2032, 3 V Lithium
Surface	Stud texture design	Type of protection	IP54
Weight of persons	At least 10 kg	Ambient temperature	0 up to +55° C
Chemical resistance	Water (room temp.), body fluids, customary disinfectants, max. 70 vol% alcohol	<i>Subject to technical modifications</i>	

9 EU Declaration of Conformity

See attachment

10 WEEE

Devices with this symbol must be treated separately during disposal. This must be done in accordance with the laws of the respective countries for environmentally sound disposal, processing and recycling of electrical and electronic equipment.

11 Contact

BBC Bircher Smart Access, BBC Bircher AG, Wiesengasse 20, CH-8222 Beringen, www.bircher.com

Designed in Switzerland / Made in Switzerland

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


CareMat A01T B434
CareMat A11T B434
1100 x 700 mm

CareMat B01T B434
CareMat B11T B434
700 x 400 mm

CareMat C01T B434
CareMat C11T B434
1100 x 700 mm,
demi-rond

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